



User Guide

Foot & Ankle for Joint Injection

Product No: 30100



MICE SWISS Sàrl
Industriestrasse 177 / 26
CH-3178 Böisingen

Tél. +41 (0) 78 727 01 44 Tél. +41 (0) 78 727 74 75

contact@mice-swiss.com



Foot & Ankle for Joint Injection

30100

This brand new model is a welcome addition to the existing Limbs & Things joint injection range of products. Like the Shoulder, Elbow and Hand & Wrist its purpose is for practising soft tissue injection for the treatment of injuries and arthritis. This product has been designed in close collaboration with Mr Jim Pickard and Dr Chris Davies from the University of Huddersfield, UK.

Skills

- Patient posture and management
- Palpation techniques
- Identifying anatomical landmarks
- Injection sites for:
 - metatarsophalangeal (MTP) joint
 - Morton's neuroma
 - tarsal tunnel
 - plantar fascia
 - sinus tarsi
 - retrocalcaneal bursa
 - tibiotalar joint
 - subtalar joint

Features



- Latex free
- Represents normal anatomy as required for palpation purposes
- Easy-to-use Feedback Console: LEDs light up when correct injection sites have been acquired

Package supplied

- **30095** 1 Foot & Ankle Needle Set
- **30101** 1 Foot & Ankle Skin
- **70102** 1 Foot & Ankle Musculature
- **BKS1846195004** Joint and Soft Tissue Injection Book
- 1 Foot & Ankle base anatomy
- 1 Feedback Console with batteries
- 1 Carry case

MICE SWISS Sàrl

Industriestrasse 177 / 26
CH-3178 Böisingen

Tél. +41 (0) 78 727 01 44 Tél. +41 (0) 78 727 74 75

contact@mice-swiss.com

